

Rehabilitation Aids

▶ Crutches

India Medico's crutches are light weight but strong. Available in wooden and Aluminium models, the crutches are quality Rehabilitation Aids appreciated and imported by several countries in the World.



IMI - 3051: Crutches Underarm (Wooden, Adjustable)



Special light weight but strong wood construction. Adjustable in height in 2.5cm increment. Fitted with foam padded under arm piece and mounted on special broad based non-slip tips. Available in three sizes, youth, adult and tall adult.

[\[Request Quotation / More info...\]](#) 

IMI - 3052: Crutches Underarm (Aluminium, Adjustable)



Lightweight P.V.C. sheathed strong aluminium pipe crutches are fitted with high density rubber underarm support & hand grip. Height & handgrip adjustments are easily made with wing nuts in one inch increments. Fitted with rubber shoes. Available in four sizes i.e. 40", 44", 48", & 52" height, with Six inch height adjustments.

P.V.C. sheathing is available in Red, Blue, Yellow or Black colours.

Packing Details:

CARTON: 54" x 18" x 18" for 20 pairs.

[Request Quotation / More info...] 

IMI - 3055: Crutches Elbow / Forearm



Telescopic, light weight aluminium tube constructed unibody crutch has fore-arm length & hand to ground distance adjustable by spring locking pins in proper notch. Forearm is fitted with plastic cuff. Fitted with rubber ferule and mounted on board base P.V.C. tips. Epoxy22 powder coated finish.

[Request Quotation / More info...] 

IMI - 3070: Tips for Walking Cane / Crutch (19 mm, 22 mm, 25 mm diameter)



Broad based non- slip moulded tip, reinforced with a steel insert that is moulded into the base. Suitable for cane and crutches. Available in 3 internal sizes of 19 mm, 22 mm & 25 mm.

[\[Request Quotation / More info...\]](#) 🏠