

## Panchakarma Products

### ▶ Dravaya Swedan Yantra (Water / Oil Bath Tables)

We have quality products in Oil Bath Tables. All our models are efficient and display quality in Panchakarma (Ayurvedic) products.



#### **IMI - 2255: FULL- BODY OIL BATH (ELEC.) AVAGHAN-SWEDAN (STAINLESS STEEL)**

Stainless Steel constructed, double walled & well insulated bath tub is 48"long x 22"wide x 20"deep. Tub is fitted with 1Kw heater and a thermostate to maintain the temprature for soothing and relaxant hydro-massage bath.



A MOTOR PUMP is provided for with a hand nozzel for specific jet massage or store out the Water or Oil. The well polished tub is fitted with patient hold bars and bottom is shaped for a proper sit of the patiient. Tub comes complete with water Inlet, Outlet and Over-flow valves.

[\[Request Quotation / More info...\]](#) 🏠

---

#### **IMI - 2501: HYDROTHERAPY TANK** (Table Shape)

Full body 'T' shape hydrotherapy tank for hydro-massage. The heavy gauge stainless steel tank, rests on 40cm (16") high tubular stand re-inforced with cross bars.

Overall size of the tank is 190cm (75") x 185cm (73") x 93cm (37").



Width of the narrow portion of the tank is 91cm (36") and depth inside the tank is 48cm (19"). Comes with Two motorized electric turbine ejectors & aerators with turbine carriages and elevators, Digital Thermometer, adjustable Head-Rest, Body Stretcher and Sitting seat fitted with nylon cloth.

Tank is fitted with Two water inlets, One overflow and one Drain fittings. Tank is finished bright inside & painted outside.

[\[Request Quotation / More info...\]](#) 🐾

---

### **IMI - 2507: WHIRLPOOL BATH** (Large Size)



Extremity Whirlpool is long, shallow & accessible to fit a variety of individual needs. Designed to treat the entire body, the unit can be beneficial for those suffering from Hip & Back conditions as well as upper Back, Neck & Shoulder problems.

#### **SPECIFICATIONS:**

-Construction: 1.6mm thick Stainless Steel Sheet Tank is continuously welded with

Argon Gas Arc welding & mounted on 10cm dia. 4 heavy duty casters.

-Turbine: Turbine Ejector-Aerator is powered by a vertically mounted 1/2 hp elec. motor & mounted on a spring loaded turbine elevator for height adjust.

-Size: 120cm x 60cm x 60cm inside the Tank.

-Water Heater: 2kw. special Heater, with cover.

-Temp. Control: Digital Temperature Controller cum Thermometer.

-Drainage: Drainage system for emptying.

-Accessories: Anodized Aluminium adjustable Seat with Arm rest.

-Power Supply: 220V AC.

[\[Request Quotation / More info...\]](#) 🏠