

Panchakarma Products

▶ General Treatment Equipment

We have quality products in General Treatment Equipment. All our models are efficient and display quality in Panchakarma (Ayurvedic) products.



IMI - 3111: EXAMINATION COUCH (Wooden)



Foam rubber padded top, size 72" x 24", upholstered with rexine & had adjustable back rest of 18" length which can be adjusted in height for desired position.

Two drawers and one writing board is provided underneath the top. Top is mounted on strudily built well polished hardwood legs.

Total height 31".

[\[Request Quotation / More info...\]](#)

IMI - 3115: HIGH-LOW TREATMENT TABLE (Electrical)



Electric, variable height, three section table is suitable for all treatments, including manipulation Both ends elevats so that the back or legs can be raised without turning the patient. One end also lower 20 deg. for some postural drainage & manipulation techniques and more comfortable sitting.



The foam padded top is 195 cm long and 70 cm wide. Table height is adjustable from 50 cm to 90 cm.

Frame-work in light blue oven baked finished.

[\[Request Quotation / More info...\]](#)

IMI - 3120: TILT TABLE (Manual)

Manually operate Tilt Table comes with Three wide straps to hold the patient. Works to standing position manually with worm gear arrangement. Range of tilt is calibrated from 0-90 deg.



The upholstered table top is 61cm wide x 195cm long x 80cm high, and is fitted on heavy duty 35mm square steel tube legs. Mounted on four lockable 10cm dia. wheels for mobility. Table comes complete with an Utility Tray & Two Gripping Handles for A.D.L. activities.

Frame beautifully finished in oven baked light blue shade.

[\[Request Quotation / More info...\]](#) 

IMI - 3141: VARUN CONTINUOUS PASSIVE MOTION UNIT

(Large, Electrical)

The CPM is designed for lower extremities (Knee & Hip), to stimulate the healing of cartilage in addition to the reduction of post operative swelling and pain. The solid state circuitry control panel with highly visual LED Digital readouts allows Flexing and Extending a joint through a pre-fixed range of motion and time. Proper fit of Calf and thigh lengths are easily achieved. The necessary limb support provides a comfortable surface for prolonged contact and are washable in cold water. The light weight unit is constructed of Aluminium Sheet and is fitted with carrying Handle and Bed fixing clamps.



SPECIFICATIONS:

-Flexion Angle: Digital display of Flexion & Extension Angles. 0 to 108 deg.

-Timer: Digital Timer upto 99 minuts.

-Angle Adjustment: Two Inching knobs to adjust starting angle.

-Speed: Fixed single speed. Right or Left Leg can be exercised by little adjustments.

-Patient Safety: Patient safety switch provided.

[\[Request Quotation / More info...\]](#) 

IMI - 3142: PANKAJ CONTINUOUS PASSIVE MOTION UNIT (U/E)



The PANKAJ CPM is designed for upper extremities (Elbow), to stimulate the healing of cartilage in addition to the reduction of post operative swelling & pain. The solid state circuitry control panel with highly visual LED digital readouts, allows the Flexing & Extending an Elbow joint through a pre-fixed range of Motion & Time. Wrist rotation and Wrist roll movement can also be achieved. Unit is fitted on mobile Electrical Jack to adjust height according to patient's requirement. The necessary limb supports provide a comfortable surface for prolonged contact. Comes with necessary accessories for Supination & Pronation, Wrist Roll and Elbow passive motions.

SPECIFICATIONS:

-Flexion Angles: Adjustable from 0 - 160 degree & displayed by Bar-graph.

-Inching: Flexion or Extension Starting position can be set by switches. -Timer: Digitally displayed, 0 - 99 minuts.

-Height Adjustment: Electrically, by Push button switches. Right or Left

arm can be exercised by little adjustments.

Safety: Patient hand held safety switch.

[\[Request Quotation / More info...\]](#) 

IMI - 3153: VIBRATOR / MASSAGER (Electrical)

No picture available

Hand-held portable, heavy duty, vibrator with wide base.

Works on 220V AC.

[\[Request Quotation / More info...\]](#) 

IMI - 2592: ULTRASOUND THERAPY UNIT (Solid State)

The I.C. based high power IMI-SONIC ultrasound therapy unit gives you a choice of treatment modes Continuous or in Pulsed wave forms. The compact IMI-SONIC with built in of 0-15 minuts gives a wide range of treatment situations. The U.S output can be adjusted by intensity knob and can be read from the output meter in watts per square centimeter.



SPECIFICATION:

-OUTPUT POWER: 15 watts continuously variable : 21 watts pulsed.

-OPERATING FREQUENCY: 1 MHz.

-TREATMENT HEAD: One.

-EFFECTIVE AREA OF RADIATION: 6 square centimeters.

-TREATMENT TIMER: 0 to 15 minutes.

-ACCESSORIES: Ultra-sound Gel

bottel.

-POWER REQUIRED: 220V.AC, 50 Hz, 100 watts.

[\[Request Quotation / More info...\]](#) 

IMI - 3155: COMPRESSION THERAPY UNIT (MEDOMER)



Medomer air massager simulates the action of muscular contraction and relaxation by repetitive inflation and deflation of the air chambers in automatic controlled sequence thereby promoting lymphatic and venous circulation.



The unit works on A.C.220V, and a linear air compressor features low noise and easy maintenance. Pressure range is adjustable between 60-100 mmHg. An automatic timer is fitted on the panel to prefix the duration of patient treatment.

The unit comes complete with Air Compressor, One Lower Limb Cuff, One Upper Limb Cuff and Two Extension Zippers.

[\[Request Quotation / More info...\]](#) 