

Physiotherapy Equipments

▶ Hydro Therapy Equipments

We have quality products in Hydro Therapy Equipments. All our models are efficient and display quality in Physiotherapy Equipments.



IMI - 2501: HYDROTHERAPY TANK (T. Shape)



Full body 'T' shape hydrotherapy tank for hydro-massage. The heavy gauge stainless steel tank, rests on 40cm (16") high tubular stand re-inforced with cross bars.

Overall size of the tank is 190cm(75") x 185cm (73") x 93cm (37").

[\[Request Quotation / More info...\]](#) 📄

Width of the narrow portion of the tank is 91cm (36") and depth inside the tank is 48cm (19"). Comes with Two motorized electric turbine ejectors & aerators with turbine carriages and elevators, Digital Thermometer, adjustable Head-Rest, Body Stretcher and Sitting seat fitted with nylon cloth. Tank is fitted with Two water inlets, One overflow and one Drain fittings.

Tank is finished bright inside & painted outside.

IMI - 2502: HYDROTHERAPY TANK

pipes and all necessary piping and fitting to wall.

IMI - 2503: ACCESSORIES FOR HYDROTHERAPY TANK



The optional accessories required for Hydrotherapy tank for conveniently handling of patients.

1. Hoist:

An Electric Hoist of 500 kg. capacity for lifting and lowering of patient in the tank. Works on 440 V. three phase A.C. power supply. Hoist is fitted with electric Travelling Trolley that moves the hoist side-ways on the 'I' beam. Comes complete with S.S cross bar with eye hooks for stretcher.

2. Trolley:

Painted trolley frame, mounted on four 20cm. dia wheels for carrying stainless steel Stretcher with nylon clothe & belts.

[\[Request Quotation / More info...\]](#) 

NOTE:

'I' beam of 8"x 4" size will provided by the department and got fixed 13 to 15 feet high from floor. Length of beam depends upon the size of the room. Minimum working length 12 feet.

IMI - 2506: WHIRLPOOL BATH (for Arm, Foot & Lower Leg)



Heavy gauge stainless steel tank sturdily fabricated and bright finished. Provides effective hydro-massage of arms, hips, legs, knees and feet. Mounted on heavy duty 10cm dia four ball bearing castors. Comes complete with removable inside aluminium Seat & Arm rest. Unit is fitted with a motorized turbine ejector and aerator on a spring loaded turbine elevator, digital thermometer, thermostate & immersion heater. Provided with drain pipe for emptying.

Inside measurements of tank 90cm long x 50cm wide x 70cm deep.
No plumbing required.

[\[Request Quotation / More info...\]](#) 

IMI - 2507: WHIRLPOOL BATH (Large Size)



Extremity Whirl-pool is long, shallow & accessible to fit a variety of individual needs. Designed to treat the entire

