

## Panchakarma Products

### ▶ Snehan Karma

We have quality products in Snehan Karma. All our models are efficient and display quality in Panchakarma (Ayurvedic) products.



#### **IMI - 2248: MASSAGE cum SHIRODHARA TABLE**



Hardwood constructed with natural wood polish finish, massage table has a slanted and curved Teak-wood top for full body oil message. Table has a shelf for storage of Oil bottles and Towels etc. Head end is fitted with a stainless steel sink for oil collecting and an overhead rod with oil pot is also provided for Shirodhara.

Approximate overall size of the table: 90" x 32"x 32" high.

#### **Packing Details:**

Crate size: 95" x 37" x 11" high - 1.

Carton size: 21" x 20" x 12" high - 1.

[\[Request Quotation / More info...\]](#) ^

#### **IMI - 2271: MASSAGE TABLE with SHIRODHARA YANTRA**

(Electrical)



Polished Teak-wood constructed table is 2 meter long and 75cm wide, has a little slant towards foot end. Yantra is fitted with stainless steel sink for Oil Collecting, Filtration and Thermostatic Heat Control arrangements of the oil. Oil is lifted and Tel Dhara drops on patients forehead automatically. Flow of Oil Drops and Movement of oil pipe can be controlled electrically. An electrical Timer Stops the treatment at the pre-fixed Time.

[Request Quotation / More info...] 

---

**IMI - 2272: SARVANG DHARA SNEHAN YANTRA** (Electrical)



Specifically designed, deluxe model Pizhichil, yantra provides Full body oil bath & Shirodhara oil bath with individual electrical controls.

**SPECIFICATION:** Construction.

**\* Dhara Patti. :**

Wooden table top with Stainless (Table) steel sink is made of One piece Imported board having 5cm high rim all around & mount on carved hard-wood legs.

**\* Dhara System.**

1. Sarvang dhara frame is made of stainless steel tube with holes for oil flow on the body.
2. Shirodhara system, made of S.S. tube is adjustable in height and swings horizontally on the forehead automatically.

**\* Snehan yantra. :**

Yantra is fitted with, low voltage D.C.motors for pumps & Systems, Oil filtration and oil Collecting chamber is fitted with Thermostatic temprature controle arrangements of the oil.

**SPECIAL FEATURES:-** The oil is automatically lifted and drops on the forehead or on full body. Flow of the oil may be controlled digitaly. An electrical Timer stops the treatment at pre-fixed time. A heater keeps the oil soothing hot for the treatment at the pre-fixed temperature. The yantra is elegantly finished.

[\[Request Quotation / More info...\]](#) 

**IMI - 2286: HI-LOW MESSAGE cum SHIRODHARA TABLE with Shirodhara Yantra** (Electrical)



Variable height dharapathy is suitable for Abhyanga, Shirodhara and Nasyam etc. The top is electrically adjustable in height with imported jack for Sit-to-Stand treatments & an easy access & transfer of patients from wheel chairs. The top is mounted on metal Frame-work with light blue / brown oven baked finish.

**FEATURES:-**

**\* Height Adjustment. :**

Electrically, with Hand- switch, from 52cm to 95cm.

**\* Dharapathy/Abhyanga :**

Size 210cm long x 80cm wide -pattika. Constructed of hard wood board with natural finish.

**\* Shirodhara Yantra :**

Yantra is fitted with Oil Collecting, Filtration & thermostatic Heat Control arrangements of the oil. Oil is lifted & Tel Dhara drops on patients forehead automatically. Flow of Oil Drops and Movement of oil nozzle can be controlled electrically. An electrical Timer Stop the treatment at the pre-fixed Time.

**\* Power Supply. : 220V AC.**

[\[Request Quotation / More info...\]](#) 

---

**IMI - 2288: HI-LOW MASSAGE cum SHIRODHARA TABLE with Shirodhara Stand (Manual)**



Variable height dharapathy is suitable for Abhyanga Shirodhara and Nasyam etc. The top is electrically adjustable in height for sit-to-stand treatments, & an easy access & transfer of patients from wheel- chairs. The top is mounted on metal Frame-work with light blue / brown oven baked finished.

**FEATURES:-**

**\* Height Adjustment. :**

Electrically, with Hand- switch, from 52cm to 95cm.

**\* Dharapathy/Abhyanga :**

Size 210cm long x 80cm wide -pattika. Constructed of hard wood board with natural finish.

**\* Shirodhara Stand :**

Painted metal stand with two wheels & S.S oil pot.

**\* Power Supply. :** 220V AC.

[\[Request Quotation / More info...\]](#) 

---

**IMI - 2283: STAND for SHIRODHARA**



Painted metal stand with chrome plated over-head bar is mounted on a wide tubular base, base is fitted with two lockable wheels for mobility. Unit comes with Stainless steel oil-pot & hanging chains.

[\[Request Quotation / More info...\]](#) 

---

**IMI - 2284: STAND for SHIRODHARA (Wooden)**



Heavy duty carved Wooden stand with wide wooden base for stability. Base is mounted on four wheels for mobility. Unit comes with "Copper" oil pot and hanging chains.

[\[Request Quotation / More info...\]](#) 

---

**IMI - 2285: SHIRODHARA YANTRA (Elec.)**



Shirodhara Yantra is fitted with, low voltage D.C.motors for Pump & System, Oil filtration and Oil Collecting Stainless steel chamber is fitted with Thermostatic temperature control arrangements of the oil.

**SPECIAL FEATURES:**

- The oil is automatically lifted and drops on patient's forehead.
- Flow of the oil may be controlled digitally.
- Movement of oil nozzle can be controlled electrically.
- An electrical Timer stops the treatment at pre-fixed time.
- A heater keeps the oil soothing hot for the treatment at the pre-fixed temperature.
- The yantra is mounted on wheels for mobility.

[\[Request Quotation / More info...\]](#) 

---