

Rehabilitation Aids

▶ Walkers / Rollator

We have widest range of Walkers / Rollators suitable for adults, children and infants. Light weight, comfortable and mostly oven baked finished. All our models are efficient and display quality in Walking Aids.



Walker Invalid's



Tubular steel frame mounted on 100 m.m. size ball bearing rubber castors. Height of under arm supports are adjustable from 105 c.m. to 135 c.m. Comes with detachable foam padded seat and rexine made back, for rest. Oven baked finish.

[\[Request Quotation / More info...\]](#) ↗

IMI - 3021: Walker Invalid's (Adult)

IMI - 3022: Walker Invalid's (Child)

IMI - 3023: Walker Infant's



For age group of 1½ to 5 years. Chrome plated rectangular shaped ring is fitted with adjustable saddle gear to hold the child in an upright position. Height of the ring and the handle provided at the rear for the patient to hold, is adjustable. Mounted on four 50mm size ball bearing rubber castors. Oven baked finish.

[\[Request Quotation / More info...\]](#) ↗

IMI - 3026: Prone Crawler



Aids elementary crawling when body support is required. Rexine made body hammock tilts forward or backward for correct positioning of the patient. Height is adjustable to accommodate children from 2 to 7 years of age. Finished oven baked.

[\[Request Quotation / More info...\]](#)

Walker Scissors Gait Preventive



Tubular steel construction. Chrome plated balance ring and handle bars are adjustable in height. Entrances from rear through hinged balance rings, height of frame is 70 cm, height of balance ring adjustable from 80 cm to 100 cm or as desired. Fitted with a bent pipe below to separate legs for scissors gait. Mounted on four 75 mm size ball bearing rubber castors. Oven baked finish.

[\[Request Quotation / More info...\]](#)

IMI - 3030: Walker Scissors Gait Preventive (Junior)

IMI - 3031: Walker Scissors Gait Preventive (Tiny Tot)

Walker Rollator



A light weight walker, eliminate side motion. Tubular steel construction with two 10 cm size fixed wheels in the front and rubber tips on the rear legs. Each handle bar is individually adjustable in height for an exact fit. By pushing with slight pressure on the handle bar, patient can move forward with safety. Oven baked finish. Child size, adjusts from 52cm to 65cm.

[\[Request Quotation / More info...\]](#)

IMI - 3033: Walker Rollator (Child)

IMI - 3035: Walker Rollator (Adult)

IMI - 3038: Walker Paraplegia (Adult)

Strong tubular construction. This walking aid provides complete balance support to a paraplegia patient and is fitted with shock absorbers. Comes with following replaceable and adjustable accessories to suit various patient.

[a] Pelvic harness with a belt and fastner to hold the patient firmly.

[b] Replaceable under arm crutches.

[c] Replaceable handle in front to hold with both hand.

[d] Replaceable can type handles.

[e] Replaceable forearm rest with hand grip.

[f] Replaceable seat for rest.

[\[Request Quotation / More info...\]](#) ↗

IMI - 3042: Walking Aid (Aluminium, Fixed Frame & Height)

Light weight anodized, polished aluminium tubular construction. The fixed frame walker is fitted with soft hand grips and mounted on P.V.C. broad based tips. Standard fixed height is 90 c.m.

[\[Request Quotation / More info...\]](#) ↗

Walking Aid Folding & Adjustable Height



Light weight, painted aluminium tubular construction, folds easily and quickly. Height of frame is adjustable from 75 c.m. to 90 c.m. by slipping locking pins in proper notch. Fitted with soft PVC hand grips and shoes.

Packing Details:

Single piece in a CARTON & Five pieces in a master BOX.

[\[Request Quotation / More info...\]](#) ↗

IMI - 3043: Walking Aid (Aluminium, Folding Frame & Adjustable Height)

IMI - 3044: Walking Aid (Aluminium, Folding Frame with two wheels & Adjustable Height)