

Rehabilitation Aids

▶ Wheelchairs

If you are looking for high quality yet economical wheelchairs, you have come to the right place...Trust India Medico Instruments for comfortable, easy to use wheelchairs which are trusted since long and known for their quality all over the world.



IMI - 3005: Wheel Chair Folding (Adult, Painted)



Available in different color combinations.

- Frame - Tubular frame, Powder Coated finish, White.
- Seat - 17" (42cm) wide, leatherette made Seat & Back (with three supporting straps under the Seat).
- Armrest - Fixed, with Skirt Guard & Arm-supports.
- Foot-rest - Fixed, with Plastic Foot Plates.
- Front Wheels - 4" (10cm) dia. Casters.
- Rear Wheels - 24" (60cm) dia., with Solid Tyres & Rings.
- Breaks - Forward Push safety Breaks.
- Packing - Carton Packing, Size 41"x 38"x 13", One piece.

[\[Request Quotation / More info...\]](#) ↗

IMI - 3011: Wheel Chair Folding (Adult Deluxe Model)



- All bright chrome plated finish.
- Folds Easily. Runs Fast & Freely.
- Fitted with forward push safety brakes.
- Seat and back made of leatherette reinforced with heavy canvas and plastic strips.
- Detachable arm rests available in full length or desk length.
- Self adjusting swivel front wheels size 20cm.
- Rear wheels fitted with hard rubber tyres size 60cm and self propelling rings.
- Autolocking, Swinging and Detachable broad foot rests. Fitted with leg support band. Heel support can be fitted on foot rests if orderd.

[\[Request Quotation / More info...\]](#) ↗

IMI - 3014: Wheel Chair Folding (With Commode)



Adult, folding wheel chair with commode seat is fitted with toggle style brakes for bedside & over the toilet use. The chair frame work is sturdily constructed from ERW steel tubing with oven baked finish. Chair is mounted on two rear wheels of 60cm dia. and two front wheels of 20cm dia., and wheels are fitted with hard rubber tyres. Foam padded Arms easily detached for entering and leaving the chair from either side. Chair comes complete with padded seat cover and commode pan.

[\[Request Quotation / More info...\]](#) ↗